

Design and Synthesis of Pyrrolidinyl Ferrocene-containing Ligands and Their Application in Highly Enantioselective Rhodium-catalyzed Olefin Hydrogenation










Xin Li, Cian Kingston, Yannick Ortin and Patrick J. Guiry*









Supplementary Information

Table of Contents


Methods for the Determination of Enantiomeric Excess (Table S1)	2
SFC Chromatograms of Racemic and Enantioenriched Compounds.....	4
¹ H NMR, ¹³ C NMR and ³¹ P NMR Spectra of New Compounds.....	21
X-Ray Crystallographic Structure for L5	49

Methods for the Determination of Enantiomeric Excess (Table S1)



Entry	Product	Assay Conditions	Retention time (min) of major enantiomer	Retention time (min) of minor enantiomer
1	 17a	HPLC, Chiralcel AD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	11.59	9.07
2	 17b	HPLC, Chiralcel AD column, Isocratic 85/15 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	10.47	8.14
3	 17c	HPLC, Chiralcel AD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	11.68	9.07
4	 17d	HPLC, Chiralcel AD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	12.44	10.47
5	 17e	HPLC, Chiralcel AD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	12.58	10.88
6	 17f	HPLC, Chiralcel AD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 254nm	29.97	27.33
7	 17g	HPLC, Chiralcel AD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	10.42	8.89
8	 17h	HPLC, Chiralcel OD column, Isocratic 95/5 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	24.03	21.02
9	 17i	HPLC, Chiralcel AD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	12.58	9.41



10	 <p>17j</p>	HPLC, Chiralcel OD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	13.25	10.53
11	 <p>17k</p>	HPLC, Chiralcel AD column, Isocratic 93/7 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 210nm	10.83	9.45
12	 <p>19a</p>	HPLC, Chiralcel AD column, Isocratic 94/6 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	11.63	9.72
13	 <p>19b</p>	HPLC, Chiralcel AD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 254nm	9.95	7.63
14	 <p>19c</p>	HPLC, Chiralcel AD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 254nm	11.13	9.28
15	 <p>19d</p>	HPLC, Chiralcel AD column, Isocratic 93/7 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	15.18	12.57
16	 <p>19e</p>	HPLC, Chiralcel OD column, Isocratic 90/10 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 210nm	9.50	12.83
17	 <p>19f</p>	HPLC, Chiralcel OD column, Isocratic 95/5 IPA/ <i>n</i> -hex, 25 °C oven, 1 mL/min flow, 220nm	30.91	25.78

SFC Chromatograms of Racemic and Enantioenriched Compounds




17a







Detector A (220nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	8.142	16695674	48.577	747180	55.722
2	10.117	17673937	51.423	593719	44.278
Totals		34369611	100.000	1340899	100.000



Detector A (220nm)


Pk #	Retention Time	Area	Area %	Height	Height %
1	10.467	46948277	100.000	1099114	100.000
Totals		46948277	100.000	1099114	100.000



Detector A (220nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	8.617	3224503	46.140	132105	56.177
2	11.175	3764000	53.860	103054	43.823



Totals		6988503	100.000	235159	100.000
--------	--	---------	---------	--------	---------



Detector A (220nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	9.067	168652	0.299	5646	0.461
2	11.683	56254004	99.701	1218756	99.539


Totals		56422656	100.000	1224402	100.000
--------	--	----------	---------	---------	---------



Detector A (220nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	10.792	59289906	47.913	1798402	52.546
2	12.892	64456210	52.087	1624154	47.454


Totals		123746116	100.000	3422556	100.000
--------	--	-----------	---------	---------	---------




Detector A (220nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	10.467	221577	0.285	6516	0.392
2	12.442	77600571	99.715	1655071	99.608

Totals		77822148	100.000	1661587	100.000
--------	--	----------	---------	---------	---------




17e





Detector A (220nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	11.192	5755130	48.120	206315	55.196
2	14.233	6204842	51.880	167472	44.804
Totals		11959972	100.000	373787	100.000




Detector A (220nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	9.975	63140	0.395	2903	0.576
2	10.875	122709	0.768 ✓	3681	0.730
3	12.575	15797125	98.837 ✓	497357	98.693
Totals		15982974	100.000	503941	100.000




Detector A (254nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	28.567	16899144	49.918	251033	52.164
2	32.217	16954378	50.082	230205	47.836
Totals		33853522	100.000	481238	100.000




Detector A (254nm)


Pk #	Retention Time	Area	Area %	Height	Height %
1	27.325	499490	0.716	7756	0.925
2	29.967	69265058	99.284	830569	99.075
Totals		69764548	100.000	838325	100.000




17g




Detector A (220nm)					
Pk #	Retention Time	Area	Area %	Height	Height %
1	8.958	21923789	48.043	861946	52.980
2	10.533	23710224	51.957	764995	47.020
Totals		45634013	100.000	1626941	100.000




Detector A (220nm)					
Pk #	Retention Time	Area	Area %	Height	Height %
1	8.892	138584	0.276	4925	0.395
2	10.417	49997588	99.724	1242624	99.605
Totals		50136172	100.000	1247549	100.000




17h




Detector A (220nm)					
Pk #	Retention Time	Area	Area %	Height	Height %
1	23.317	12842045	48.218	227695	50.932
2	29.783	13791122	51.782	219366	49.068
Totals		26633167	100.000	447061	100.000




Detector A (220nm)					
Pk #	Retention Time	Area	Area %	Height	Height %
1	✓ 21.017	726459	0.486 ✓	17448	0.873
2	21.833	717214	0.480	15123	0.756
3	24.033	148073711	99.034 ✓	1996872	98.371
Totals		149517384	100.000	1999443	100.000




17i




Detector A (220nm)					
Pk #	Retention Time	Area	Area %	Height	Height %
1	9.533	36998543	49.678	1346372	57.003
2	12.750	37478706	50.322	1015557	42.997
Totals		74477249	100.000	2361929	100.000



Detector A (220nm)					
Pk #	Retention Time	Area	Area %	Height	Height %
1	9.408	130207	0.159	4008	0.240
2	12.583	81890223	99.841	1664224	99.760
Totals		82020430	100.000	1668232	100.000




17j






Detector A (220nm)


Pk #	Retention Time	Area	Area %	Height	Height %
1	14.258	17461334	49.827	505483	56.831
2	19.200	17582931	50.173	383966	43.169
Totals		35044265	100.000	889449	100.000




Detector A (254nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	10.525	16279	0.423	577	0.802
2	13.250	3827775	99.577	71342	99.198
Totals		3844054	100.000	71919	100.000







19a





Detector A (220nm)					
Pk #	Retention Time	Area	Area %	Height	Height %
1	9.833	24847348	49.782	759934	52.739
2	12.225	25065435	50.218	681002	47.261
Totals		49912783	100.000	1440936	100.000




Detector A (220nm)					
Pk #	Retention Time	Area	Area %	Height	Height %
1	9.717	565978	1.791	18055	2.672
2	11.633	31043172	98.209	657640	97.328
Totals		31609150	100.000	675695	100.000




19b






19c



Detector A (254nm)

Pk #	Retention Time	Area	Area %	Height	Height %
1	8.667	13432742	49.586	375088	56.105
2	11.517	13657216	50.414	293462	43.895




Totals		27089958	100.000	668550	100.000
--------	--	----------	---------	--------	---------




Detector A (220nm)



Pk #	Retention Time	Area	Area %	Height	Height %
1	9.275	866539	1.186	24482	1.666
2	11.125	72205768	98.814	1444973	98.334


Totals		73072307	100.000	1469455	100.000
--------	--	----------	---------	---------	---------







19e








19f




¹H NMR, ¹³C NMR and ³¹P NMR Spectra of New Compounds





¹H NMR Spectrum of (*R,S_p*)-14




^{13}C NMR Spectrum of (*R,Sp*)-**14**




^{31}P NMR Spectrum of (*R,Sp*)-**14**





¹H NMR Spectrum of (*R,R*)-**14**




^{13}C NMR Spectrum of (*R,Rp*)-**14**




^{31}P NMR Spectrum of (*R,Rp*)-**14**





¹H NMR Spectrum of (*R,S_p*)-15




^{13}C NMR Spectrum of (*R,S*)-**15**




^{31}P NMR Spectrum of (*R,S*)-**15**





¹H NMR Spectrum of (*R,R_p*)-15




^{13}C NMR Spectrum of (*R,Rp*)-**15**




^{31}P NMR Spectrum of (*R,Rp*)-**15**




¹H NMR Spectrum of ligand L1




^{13}C NMR Spectrum of ligand **L1**




^{31}P NMR Spectrum of ligand **L1**





^{31}P - ^{31}P COSY of Ligand **L1**




L2




¹H NMR Spectrum of ligand L2







L3




¹H NMR Spectrum of ligand L3




^{13}C NMR Spectrum of ligand **L3**




^{31}P NMR Spectrum of ligand **L3**





^{19}F NMR Spectrum of ligand **L3**




L4





¹H NMR Spectrum of ligand L4






^{19}F NMR Spectrum of ligand **L4**







^{13}C NMR Spectrum of ligand **L5**




^{31}P NMR Spectrum of ligand **L5**








^{13}C NMR Spectrum of ligand **L6**




^{31}P NMR Spectrum of ligand **L6**




¹H NMR Spectrum of ligand L7






^{13}C NMR Spectrum of ligand **L7**



^{31}P NMR Spectrum of ligand **L7**








^{31}P - ^{31}P COSY of ligand L8

X-Ray Crystallographic Structure of L5



L5 (molecule; thermal ellipsoids are drawn on the 50% probability level; O1 15% occupied)

Crystal data and structure refinement for **L5**.

Identification code	gui53
Empirical formula	C ₄₆ H ₃₇ Fe N O _{2.15} P ₂
Formula weight	755.91
Temperature	100(2) K
Wavelength	1.54184 Å
Crystal system	Monoclinic
Space group	P 21
Unit cell dimensions	a = 8.0501(1) Å b = 15.2402(3) Å c = 14.8419(2) Å
Volume	1806.83(5) Å ³
Z	2
Density (calculated)	1.389 Mg/m ³
Absorption coefficient	4.510 mm ⁻¹
F(000)	786.4
Crystal size	0.3073 x 0.2633 x 0.0342 mm ³
Theta range for data collection	4.17 to 65.19°.
Index ranges	-9 ≤ h ≤ 9, -17 ≤ k ≤ 17, -17 ≤ l ≤ 17
Reflections collected	14880
Independent reflections	6156 [R(int) = 0.0345]
Completeness to theta = 65.19°	99.8 %
Absorption correction	Analytical
Max. and min. transmission	0.860 and 0.397
Refinement method	Full-matrix least-squares on F ²
Data / restraints / parameters	6156 / 1 / 474
Goodness-of-fit on F ²	1.066
Final R indices [I > 2σ(I)]	R1 = 0.0355, wR2 = 0.0891
R indices (all data)	R1 = 0.0386, wR2 = 0.0906
Absolute structure parameter	-0.017(3)
Largest diff. peak and hole	0.247 and -0.340 e.Å ⁻³